

## Newest specialty HCD — xm18s | High Shear...



- ◆ Xennovate changes rules of game, again.
- ◆ Name chosen specifically to position our company as the innovator in hydrocolloids
- ◆ Our thin integrated formulation enables nearly cost competitive skin attachment devices that are highly skin friendly
- ◆ Nearby image shows test method for shear
- ◆ xm18s stays in place for over 20 minutes
- ◆ “Standard” HCDs, including our 18w and 18o remain in place for only about 15 seconds
- ◆ By this method, xm18s has amazingly over 80x shear adhesion of “standard” HCDs
- ◆ Interestingly, xm18e has six times standard adhesion under this method

